

neosanit® pre-soaped sponges are a cosmetic product for daily personal hygiene. They are especially indicated for a correct cleansing of the skin.

Units per package	Description	Content per box	Units per box	Measures	Sponge weight
1	Polyurethane foam of 15 kg/m <sup>3</sup>	12 bags with 100 packages	1,200 sponges	13 x 9 cm ± 5% Thickness: 2.45 cm	4,0 – 5,0 g

#### TRANSPORT INFORMATION:

Cardboard box measures: 780 x 590 x 900 mm  
 Cardboard box weight: 8 kg ± 10 %  
 Pallet with 5 cardboard boxes  
 Pallet measures: 126 x 80 x 257 cm  
 Full truck with 33 pallets

#### BAR CODES:

GTIN 13: 8426696003835  
 GTIN 14: 18426696003832

#### STORAGE AND SELF-LIFE

- Store in a cool, dry place out of direct sunlight.
- Expiry date 36 months. Once opened consume within 12 months.

#### PHYSICAL PROPERTIES

The sponges are made of polyurethane foam which provides a soft, effective, durable and lightweight material.

#### GEL: DESCRIPTION, INGREDIENTS AND PROPERTIES

Pre-soaped sponges are impregnated with dermatologic and hypoallergenic pH 5.5 ± 0.3 dye free and scent free soap.



INGREDIENTS	FUNCTION
Aqua	Solvent
Sodium Laureth Sulfate	Cleansing, emulsifying, foaming, surfactant
Cocamide DEA	Emulsifying, emulsion stabilising, foam boosting, surfactant, viscosity controlling
Disodium Cocoamphodiacetate	Cleansing, foam boosting, hair and skin conditioning, hydrotope, surfactant
Citric Acid	Buffering, chelating, masking
Lactic Acid, Propylene Glycol, Malic Acid, Rosmarinus Officinalis Extract, Camellia Sinensis Leaf Extract, Citrus Aurantium Amara Flower Extract	Preservative
Sodium Hydroxide	Buffering

## MICROBIOLOGICAL QUALITY

Limit (Aerobic bacteria, molds, yeasts) cm <sup>2</sup> of product	< 100 UFC / cm <sup>2</sup>
Escherichia coli, Staphylococcus aureus, Pseudomonas aeruginosa, Candida albicans	Absence



## INFORMATION

Neosanit® pre-soaped sponges are a cosmetic product as defined in section 201(i) of the Federal Food, Drug, and Cosmetic Act (DF&C Act) enforced by the Food and Drug Administration (FDA). Label terms are in accordance to FD&C Act and Fair Packing and Labelling Act (FP&L Act).

The product complies with California Proposition 65 January 03, 2020

The production plant is certified with UNE-EN ISO 22716:2007 regarding product documentation and Good Manufacturing Practices for cosmetic products.

The company is certified according to the following standards: UNE-EN ISO 9001:2015, UNE-EN ISO 14001:2015.

## PRECAUTIONS FOR USE

- Keep out of reach of children.
- Do not ingest.
- Avoid contact with eyes.
- For external use only.



1. The gel lathers with a very small amount of water



2. Rub the sponge gently over the skin



3. Rinse with water



4. Discard the sponge

## RECOMMENDATIONS FOR USE

- We recommend caretakers the use of protective gloves and sterile containers. Dry the skin with a drying cloth or towel.
- For bedridden patients or persons with reduced mobility, you can use a slightly damp cloth to remove the lather.
- To improve the patient's comfort, you can moisten the sponge with warm water up to 98,6 Fahrenheit



**Do not immerse  
In water**