neo sanit

Esponja con jabón de un solo uso para una correcta limpieza e higiene de la piel

neosanit® pre-soaped sponges are a cosmetic product for daily personal hygiene. They are especially indicated for a correct cleansing of the skin.

Units per package	Description	Content per box	Units per box	Measures	Sponge weight
10	BIODEGRADABLE Poly -lactic fiber of 150 g/m2 ± 5%	3 bags de 90 packages	2.700 sponges	12 x 20 cm ± 5%	3,5 – 4,5 g

TRANSPORT INFORMATION:

Cardboard box measures: 770 x 590 x 640 mm

Cardboard box weight: 10,6 kg ± 10 %

Pallet with 6 cardboard boxes
Pallet measures: 120 x 80 x 210 cm

Full truck with 33 pallets

BAR CODES:

GTIN 13: 8426696003811 GTIN 14: 18426696003818

STORAGE AND SELF-LIFE

- Store in a cool, dry place out of direct sunlight.
- Expiry date 36 months. Once opened consume within 12 months.

PHYSICAL PROPERTIES

The sponges are made of poly-lactic fiber which provides a soft, effective, durable and lightweight material.

GEL: DESCRIPTION, INGREDIENTS AND PROPERTIES

Pre-soaped sponges are impregnated with dermatologic and hypoallergenic pH 5.5 \pm 0.3 dye free and scent free soap.



neosanit

INGREDIENTS	FUNCTION	
Aqua	Solvent	
Sodium Laureth Sulfate	Cleansing, emulsifying, foaming, surfactant	
Cocamide DEA	Emulsifying, emulsion stabilising, foam boosting, surfactant, viscosity controlling	
Disodium Cocoamphodiacetate	Cleansing, foam boosting, hair and skin conditioning, hydrotope, surfactant	
Citric Acid	Buffering, chelating, masking	
Lactic Acid, Propylene Glycol, Malic Acid, Rosmarinus Officinalis Extract, Camellia Sinensis Leaf Extract, Citrus Aurantium Amara Flower Extract	Preservative	
Sodium Hydroxide	Buffering	

MICROBIOLOGICAL QUALITY	
Limit (Aerobic bacteria, molds, yeasts) cm ² of product	< 100 UFC / cm ²
Escherichia coli, Staphylococcus aureus, Pseudomonas aeruginosa, Candida albicans	Absence



neosanit

INFORMATION

The product complies with the requirements of Regulation (EC) number 1223/2009 of the European Parliament and the Council of 30 November 2009 on cosmetic products.

The skin compatibility of the product has been tested under dermatological control.

No animal testing.

The production plant is certified with UNE-EN ISO 22716:2007 regarding product documentation and Good Manufacturing Practices for cosmetic products.

The company is certified according to the following standards: UNE-EN ISO 9001:2015, UNE-EN ISO 14001:2015.

PRECAUTIONS FOR USE

- Keep out of reach of children.
- Do not ingest.
- Avoid contact with eyes.
- For external use only.

The gel lathers with

a very small amount of water

Rub the spong gently over the skin

Rinse with water

Discard the sponge

RECOMENDACIONES PARA SU APLICACIÓN

- We recommend caretakers the use of protective gloves and sterile containers.
- Dry the skin with a drying cloth or towel.
- For bedridden patients or persons with reduced mobility, you can use a slightly damp cloth to remove the lather.
- To improve the patient's comfort, you can moisten the sponge with warm water up to 98,6 Fahrenheit





TECHNICAL DATA SHEET Date of issue: MAY/2021 Edition 01 - Page 3 of 3

Camino de Fitena s/n 18240 Pinos Puente (Granada) España Tel. + 34 958 459 050 / Fax +34 958 452 300 www.jalsosa.com